# Thermovent<sup>®</sup> T Heat and Moisture Exchanger

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[www.smiths-medical.com](http://www.smiths-medical.com)**smiths medical**

# Portex® Thermovent® T Heat and Moisture Exchanger

## INSTRUCTIONS FOR USE

These instructions are intended for use with the following **Portex® product:**

Catalogue No: 100/570/015 *Thermovent® T Heat and*  
Internal Ref: 100/570/015CZ *Moisture Exchanger*

**These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.**

**NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.**

### Description

Thermovent® T is a sterile, single use heat and moisture exchanger which provides effective humidification for spontaneously breathing adult and paediatric patients whose upper airways are bypassed by a tracheal or tracheostomy tube.

### Indications for use

The Portex® Thermovent® T is indicated for:

- Preserving humidity and minimising damage to tracheal epithelial cells
- Reduction of heat loss through the broncho pulmonary tree
- Helping to prevent thickened secretions and changes in lung function during long term anaesthesia or ventilatory support.

### Instructions for use

- 1 The device should be connected to the 15mm termination of the tracheal or tracheostomy tube.
- 2 Where tenacious secretions or mucous plugging develop suitable therapy must be instituted.
- 3 The patient should be closely monitored for any complications resulting from blockage to the artificial airway or Thermovent® T. Should any complications develop, the device should be removed or replaced.
- 4 Change the device every 24 hours or more frequently as determined by expert clinical assessment to prevent accumulation of secretions in the device.

### Precautions

- 1 As with all airway management products, the security of all connections and patency of the system and device should be checked prior to use.
- 2 This device should only be used with spontaneously breathing patients.
- 3 The effect of the additional deadspace volume of this device should be evaluated on an individual patient basis.

- 4 The device must not be used on conjunction with heated or nebulized water humidification systems or nebulizers.

### Warnings

- 1 The device is intended for single use only and must not be re-used or cleaned. Wetting of the heat and moisture exchange elements with any cleaning solution may dangerously increase the resistance to breathing or may result in the retention of harmful residues or transmission of infectious agents.
- 2 Expert clinical judgement must be used in assessing the patient's humidification requirements. This device must be used as described in these instructions and in accordance with currently accepted medical techniques.
- 3 If a heat and moisture exchange element is displaced from the housing (for example by patient coughing), the device should be discarded and a new device used. Under no circumstances must replacement of the element be attempted.
- 4 The device may block if the patient develops copious secretions, pulmonary oedema or bleeding when the device is in use. If the device becomes blocked then it must be replaced. Under no circumstances should any attempt be made to clear or clean out the blockage as this may damage the heat and moisture exchange elements and could lead to further blockage.
- 5 Caution should be taken if this device is used on children as it contains components which may be removed and could be harmful if ingested or inhaled.

### Technical Specification

- 1 **Moisture Output**  
25mg/L H<sub>2</sub>O @ 10 breath/min  
Tidal Volume 1,000 ml
- 2 **Resistance to Flow**  
0.9 hPa (cmH<sub>2</sub>O) at gas flow rate 0.5L/sec (30L/min)  
2.0 hPa (cmH<sub>2</sub>O) at gas flow rate 1.0L/sec (60L/min)  
3.5 hPa (cmH<sub>2</sub>O) at gas flow rate 1.5L/sec (90L/min)

*Thermovent, Portex, Portex design mark and Smiths Medical design mark are trademarks Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners. © 2019 Smiths Medical. All rights reserved.*

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**en****Single use.**

**Do Not Reuse:** Medical devices require specific material characteristics to perform as intended. These characteristics have been verified for single use only. Any attempt to re-process the device for subsequent re-use may adversely affect the integrity of the device or lead to deterioration in performance.

**Sterile unless unit container is opened or damaged.**

**Destroy after single use. Do not resterilize.**

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Thermovent, Portex, designmerket Portex og Smiths Medical er varemærker som tilhører Smiths Medical. Symbolet® angir at varemærket er registrert ved U.S. Patent and Trademark Office og i visse andre land. Alle andre navn og merker som er nevnt, er varemærker, handelsnavn og servicemærker som tilhører de respektive eiere. © 2019 Smiths Medical. Med enerett.